

Product Datasheet

CD79b Antibody (orb1824990)

Description	CD79 (also designated Ig Alpha/Ig beta) is a heterodimer composed of Alpha chains, designated CD79A ...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	A recombinant fragment of human CD79b protein was used as the immunogen for the CD79b antibody.
Storage	Aliquot the CD79b antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Uniprot ID	P40259
Hazard Information	This CD79b antibody is available for research use only.
Dilution Range	ELISA: 2-4mg/ml for coating (order BSA-free format),Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.