

Product Datasheet

p19INK4d Cdkn2d Antibody (Monoclonal, DCS-100) (orb18251)

Description

p19INK4d Cdkn2d Antibody (Monoclonal, DCS-100)

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

ICC, IHC, WB

Applications
Immunogen

Recombinant human p19INK4d.

Form/Appearance

Lyophilized

Concentration

Adding 1 ml of PBS buffer will yield a concentration of 100 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 1ml of PBS buffer will yield a concentration of 100ug/ml.

Isotype

Mouse IgG1

Clonality

Monoclonal

MW


17809 MW

Uniprot ID
[Q60773](https://www.uniprot.org/entry/Q60773)
Dilution Range

Immunohistochemistry (Paraffin-embedded Section), 1-2µg/ml, Human, By Heat
Immunocytochemistry , 1µg/ml, Human, -Western blot, 0.5-1µg/ml, Human

Expiration Date

12 months from date of receipt.



WB analysis
of
p19INK4d
using anti-
p19IN...