

Product Datasheet

BCMA Antibody / CD269 / TNFRSF17 (orb1825307)

Description	The B cell maturation protein (BCMA) is a type I integral membrane protein that belongs to the tumor...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	A recombinant human protein fragment (within amino acids 1-184) was used as the immunogen for the B-cell maturation protein antibody.
Storage	Aliquot the B-cell maturation protein antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	Q02223
Hazard Information	This B-cell maturation protein antibody is available for research use only.
Dilution Range	ELISA (Order BSA-free format for coating),Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.