

Product Datasheet

PLK1 Antibody / Polo like kinase 1 (orb1825344)

Description	Plk (for polo-like kinase) encodes a serine/threonine kinase that is closely related to polo and CDC...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	A recombinant fragment of human protein (within amino acids 1-200) was used as the immunogen for the Polo like kinase 1 antibody.
Storage	Aliquot the Polo like kinase 1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Uniprot ID	P53350
Hazard Information	This Polo like kinase 1 antibody is available for research use only.
Dilution Range	ELISA (Order BSA-free format for coating),Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.