

Product Datasheet

SULT1E1 Antibody / Estrogen Sulfotransferase (orb1825548)

Description

The soluble sulfotransferases contribute

activation of procar...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

A recombinant partial protein sequence (within amino acids 100-294) from the human protein was used as the immunogen for the Estrogen Sulfotransferase antibody.

Storage

Aliquot the Estrogen Sulfotransferase antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note

For research use only

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG2, kappa

Clonality

Monoclonal

Uniprot ID

P49888

Hazard
Information

This Estrogen Sulfotransferase antibody is available for research use only.

Dilution Range

Western blot: 1-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT

Expiration Date

12 months from date of receipt.