

Product Datasheet

Collagen IV Antibody (orb1825581)

Description	Collagen Type IV is a major component of the basement membrane and plays an important role in cell a...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Purified human Collagen IV protein was used as the immunogen for the Collagen IV antibody.
Storage	Aliquot the Collagen IV antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P08572, P29400, P02462, P53420
Hazard Information	This Collagen IV antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.