

## Product Datasheet

**PCK Antibody (Monoclonal, C-11+PCK-26+CY-90+KS-1A3+M20+A53-B/A2) (orb18257)**

**Description**

PCK Antibody (Monoclonal, C-11+PCK-


**Species/Host**

Mouse

**Reactivity**

Human, Rat

**Conjugation**

Unconjugated

**Tested**

IHC, WB

**Applications**

Immunohistochemical  
analysis of paraffin...

**Immunogen**

mixture of several monoclonal  
cytokeratin clones.

**Form/Appearance**

Lyophilized

**Concentration**

Adding 1 ml of PBS buffer will yield a  
concentration of 100 µg/ml.

**Storage**

Store at -20°C for one year from date of  
receipt. After reconstitution, at 4°C for  
one month. It can also be aliquotted and  
stored frozen at -20°C for six months.  
Avoid repeated freeze-thaw cycles.

**Note**

For research use only

**Application notes**

Other applications have not been tested.  
Optimal dilutions should be determined  
by end users. . Add 1ml of PBS buffer will  
yield a concentration of 100ug/ml.

**Isotype**

Mouse IgG2a/IgG1

**Clonality**

Monoclonal

**MW**

66039 MW

**Uniprot ID**

**P04264**

**Dilution Range**

Immunohistochemistry (Paraffin-  
embedded Section), 0.1-0.5µg/ml,  
Human, rat, By Heat  
Immunohistochemistry (Frozen Section),  
0.1-0.5µg/ml, Human, rat, -Western blot,  
1-2µg/ml, Human, rat

**Expiration Date**

12 months from date of receipt.