


Product Datasheet

PCK-26 Krt8 Antibody (Monoclonal, PCK-26) (orb18258)


Description

PCK-26 Krt8 Antibody (Monoclonal, PCK-26)


Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Cytokeratin from human epidermis.
Form/Appearance	Lyophilized
Concentration	Adding 1 ml of PBS buffer will yield a concentration of 100 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 1ml of PBS buffer will yield a concentration of 100ug/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	54019 MW
Uniprot ID	Q10758
Dilution Range	Western blot, 0.5-1µg/ml, Human, mouse Immunohistochemistry (Paraffin-embedded Section), 1-2µg/ml, Human, By Heat Immunofluorescence, 1-2µg/ml, Human, -
Expiration Date	12 months from date of receipt.



WB analysis of Pan cytokeratin using ant...



IF analysis of Pan cytokeratin using ant...



IHC analysis of Pan cytokeratin using an...