

Product Datasheet

CD3 gamma Antibody (orb1825957)

Description	The T cell antigen receptor (TCR) recognizes foreign antigens and translates such recognition events...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 1-182) from the human protein was used as the immunogen for the CD3 gamma antibody.
Storage	Aliquot the CD3 gamma antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2c, kappa
Clonality	Monoclonal
Uniprot ID	P09693
Hazard Information	This CD3 gamma antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml
Expiration Date	12 months from date of receipt.