

Product Datasheet

DRAP1 Antibody / NC2 alpha (orb1826058)

Description	NC2 (negative cofactor 2) is a dimeric histone-fold complex that represses RNA polymerase II transcr...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FACS, IF
Immunogen	Recombinant human protein was used as the immunogen for the NC2 alpha antibody.
Storage	Aliquot the NC2 alpha antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2b, kappa
Clonality	Monoclonal
Uniprot ID	Q14919
Hazard Information	This NC2 alpha antibody is available for research use only.
Dilution Range	Flow cytometry: 1-2ug/million cells, Immunofluorescence: 1-2ug/ml
Expiration Date	12 months from date of receipt.