

Product Datasheet

MUC5AC Antibody (orb1826093)

Description	Mucins are a group of high molecular weight glycoproteins consisting of a mucin core protein and O-l...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant full-length human MUC5AC protein was used as the immunogen for the MUC5AC antibody.
Storage	Aliquot the MUC5 ^{°C} antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Mouse IgG1, kappa
Clonality	Recombinant
Uniprot ID	P98088
Hazard Information	This MUC5AC antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml
Expiration Date	12 months from date of receipt.