

## Product Datasheet

### p57 Antibody / p57Kip2 (orb1826102)

|                            |                                                                                                         |
|----------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Recognizes a protein of 57kDa, identified as p57Kip2. It shows no cross-reaction with p27Kip1. p57Ki... |
| <b>Species/Host</b>        | Mouse                                                                                                   |
| <b>Reactivity</b>          | Human                                                                                                   |
| <b>Conjugation</b>         | Unconjugated                                                                                            |
| <b>Tested Applications</b> | IHC-P                                                                                                   |
| <b>Immunogen</b>           | Recombinant full-length human p57Kip2 protein was used as the immunogen for the p57Kip2 antibody.       |
| <b>Storage</b>             | Aliquot the p57Kip2 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.    |
| <b>Note</b>                | For research use only                                                                                   |
| <b>Formula</b>             | 1 mg/ml in 1X PBS; BSA free, sodium azide free                                                          |
| <b>Isotype</b>             | Mouse IgG, kappa                                                                                        |
| <b>Clonality</b>           | Recombinant                                                                                             |
| <b>Uniprot ID</b>          | <b>P49918</b>                                                                                           |
| <b>Hazard Information</b>  | This p57Kip2 antibody is available for research use only.                                               |
| <b>Dilution Range</b>      | Immunohistochemistry (FFPE): 1-2ug/ml                                                                   |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                         |