

Product Datasheet

SI-CLP Antibody / CHID1 / Chitinase Domain Containing 1 (orb1826113)

Description	SI-CLP (stabilin-1 interacting chitinase-like), also known as CHID1 (chitinase domain-containing pro...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant full-length human SI-CLP protein was used as the immunogen for the CHID1 antibody.
Storage	Aliquot the SI-CLP antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Uniprot ID	Q9BWS9
Hazard Information	This SI-CLP antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.