


Product Datasheet

CHRAC1 antibody (orb182614)

Description	CHRAC1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human CHRAC1 (1-100/131 aa)
Target	CHRAC1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	14 kDa
Uniprot ID	Q9NRG0
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from M...