

Product Datasheet

PNLIP Antibody / Pancreatic lipase (orb1826203)

Description	The lipase gene family belongs to one of the most robust genetic superfamilies found in living organ...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	A recombinant partial protein (within amino acids 1-200) from the human protein was used as the immunogen for the Pancreatic lipase antibody.
Storage	Aliquot the Pancreatic lipase antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Uniprot ID	P16233
Hazard Information	This Pancreatic lipase antibody is available for research use only.
Dilution Range	ELISA (Order BSA-free format for coating),Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.