

Product Datasheet

KChIP2 Antibody / KCNIP2 (orb1826208)

Description	In the brain and heart, rapidly inactivating (A-type) voltage-gated potassium (Kv) currents control ...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein (within amino acids 1-270) from the human protein was used as the immunogen for the KCNIP2 antibody.
Storage	Aliquot the KCNIP2 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Uniprot ID	Q9NS61
Hazard Information	This KCNIP2 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.