

Product Datasheet

CD39 Antibody (orb1826213)

Description	CD39, also known as ectonucleoside triphosphate diphosphohydrolase 1 (ENP1), is an integral membrane...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant full-length human CD39 protein was used as the immunogen for the CD39 antibody.
Storage	Aliquot the CD39 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P49961
Hazard Information	This CD39 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.