


Product Datasheet

CDK1 Antibody (Monoclonal, PSTAIR) (orb18263)

Description	CDK1 Antibody (Monoclonal, PSTAIR)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic 16 amino acid oligopeptide containing the PSTAIR sequence conjugated to BSA.
Form/Appearance	Lyophilized
Concentration	Adding 1 ml of PBS buffer will yield a concentration of 100 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 1ml of PBS buffer will yield a concentration of 100ug/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	34135 MW
Uniprot ID	P39951
Dilution Range	Western blot, 0.5-1µg/ml, Human, mouse, rat
Expiration Date	12 months from date of receipt.



WB analysis of CDK1 using anti-CDK1 anti...