

Product Datasheet

GAD65 Antibody / GAD2 (orb1826305)

Description	This mAb recognizes a protein of 65kDa, which is identified as glutamic acid decarboxylase 2 (GDA2)....
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	A recombinant human GAD2/GAD65 protein fragment (within amino acids 6-99) was used as the immunogen for the GAD65 antibody.
Storage	Aliquot the GAD65 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	Q05329
Hazard Information	This GAD65 antibody is available for research use only.
Dilution Range	ELISA (Order BSA-free format for coating)
Expiration Date	12 months from date of receipt.