



## Product Datasheet

### Contactin 1 antibody (orb182660)

|                            |                                                                                   |
|----------------------------|-----------------------------------------------------------------------------------|
| <b>Description</b>         | Contactin 1 antibody                                                              |
| <b>Species/Host</b>        | Rabbit                                                                            |
| <b>Reactivity</b>          | Mouse                                                                             |
| <b>Conjugation</b>         | Unconjugated                                                                      |
| <b>Tested Applications</b> | FC                                                                                |
| <b>Immunogen</b>           | KLH conjugated synthetic peptide derived from human Contactin 1 (151-250/1018 aa) |
| <b>Target</b>              | Contactin 1                                                                       |
| <b>Preservatives</b>       | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.             |
| <b>Form/Appearance</b>     | Liquid                                                                            |
| <b>Concentration</b>       | 1mg/ml                                                                            |
| <b>Storage</b>             | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.  |
| <b>Note</b>                | For research use only                                                             |
| <b>Isotype</b>             | IgG                                                                               |
| <b>Clonality</b>           | Polyclonal                                                                        |
| <b>MW</b>                  | 111 kDa                                                                           |
| <b>Uniprot ID</b>          | <b>Q12860</b>                                                                     |
| <b>Dilution Range</b>      | Flow-Cyt=3ug/test, WB=1:500-2000, ELISA=1:5000-10000                              |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                   |



Blank  
control:  
Mouse brain.  
Primary  
Anti...



Blank  
control:  
Mouse  
kidney.  
Primary  
Ant...