

Product Datasheet

CLEC4C Antibody (orb1826928)

Description	Anti-CLEC4C antibody(12D6); IgG1 Chimeric mAb
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Target	CLEC4C
Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.
Note	For research use only
Application notes	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Isotype	Rabbit/Human Fc chimeric IgG1
Clonality	Recombinant
Uniprot ID	Q8WTT0
Dilution Range	Flow Cyt 1:100
Expiration Date	12 months from date of receipt.