

Product Datasheet

MMAE Antibody (orb1826938)

Description	Anti-MMAE antibody(8C4); Rabbit mAb
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Target	MMAE
Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.
Note	For research use only
Application notes	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Isotype	Rabbit IgG
Clonality	Recombinant
Uniprot ID	N.A.
Dilution Range	ELISA 1:5000-10000
Expiration Date	12 months from date of receipt.