

## Product Datasheet

### Human Interleukin 1 $\beta$ , IL-1 $\beta$ ELISA Kit (orb1826974)

|                          |                                                                  |
|--------------------------|------------------------------------------------------------------|
| <b>Description</b>       | Human Interleukin 1 $\beta$ , IL-1 $\beta$ ELISA Kit             |
| <b>Reactivity</b>        | Human                                                            |
| <b>Range</b>             | 2.344 pg/ml-150 pg/ml                                            |
| <b>Note</b>              | For research use only                                            |
| <b>Application notes</b> | Detection Wavelength: 450 nm                                     |
| <b>Sample Types</b>      | serum, plasma, tissue<br>homogenates, cell culture<br>supernates |
| <b>Uniprot ID</b>        | <b>P01584</b>                                                    |
| <b>Sensitivity</b>       | 0.586 pg/ml                                                      |
| <b>Expiration Date</b>   | Please enquire.                                                  |