

## Product Datasheet

### Mouse Interleukin 6,IL-6 ELISA Kit (orb1826976)

|                          |                                                               |
|--------------------------|---------------------------------------------------------------|
| <b>Description</b>       | Mouse Interleukin 6,IL-6 ELISA Kit                            |
| <b>Reactivity</b>        | Mouse                                                         |
| <b>Range</b>             | 3.9 pg/ml-250 pg/ml                                           |
| <b>Note</b>              | For research use only                                         |
| <b>Application notes</b> | Detection Wavelength: 450 nm                                  |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates,<br>cell culture supernates |
| <b>Uniprot ID</b>        | <b>P08505</b>                                                 |
| <b>Sensitivity</b>       | 0.975 pg/ml                                                   |
| <b>Expiration Date</b>   | Please enquire.                                               |