


Product Datasheet

CPNE8 antibody (orb182717)

Description	CPNE8 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human CPNE8 (1-100/564 aa)
Target	CPNE8
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	63 kDa
Uniprot ID	Q86YQ8
Expiration Date	12 months from date of receipt.



Western blot analysis of HT29 Cell Lysat...