

Product Datasheet

Tyrosine Hydroxylase Antibody (Monoclonal, TH-16) (orb18276)

Description

Tyrosine Hydroxylase Antibody (Monoclonal, TH-

Species/Host

Mouse

Reactivity

Human, Mouse, Rabbit, Rat

Conjugation

Unconjugated

Tested
Applications

IHC, WB

Immunogen

Rat tyrosine hydroxylase (TH).

Form/Appearance

Lyophilized

Concentration

Adding 1 ml of PBS buffer will yield a concentration of 100 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 1ml of PBS buffer will yield a concentration of 100ug/ml.

Isotype

Mouse IgG1

Clonality

Monoclonal

MW

55966 MW

Uniprot ID
P04177
Dilution Range

Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human, rat, rabbit, By Heat Immunohistochemistry (Frozen Section), 0.5-1µg/ml, Human, rat, rabbit, -Western blot, 0.25-0.5µg/ml, Human, mouse, rat, rabbit

Expiration Date

12 months from date of receipt.



Western
Blot
analysis of
rat tissue
lysa...