

Product Datasheet

Beta-Tubulin TUBB1 Antibody (Monoclonal, TUB 2.1) (orb18287)

Description

Beta-Tubulin TUBB1 Antibody (Monoclonal, TUB 2.1)

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ICC, IF, IHC, WB

Applications

Immunohistochemical
analysis of paraffin...

Immunogen

Tubulin from rat brain.

Form/Appearance

Lyophilized

Concentration

Adding 1 ml of PBS buffer will yield a concentration of 100 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 1ml of PBS buffer will yield a concentration of 100ug/ml.

Isotype

Mouse IgG1

Clonality

Monoclonal

MW

50327 MW

Uniprot ID
[Q9H4B7](#)
Dilution Range

Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human, mouse, rat, By Heat Western blot, 2µg/ml, Human, mouse, ratImmunocytochemistry/Immunofluorescence,5µg/ml, Human

Expiration Date

12 months from date of receipt.