



Product Datasheet

Cytochrome P450 2B1 antibody (orb182872)


Description	Cytochrome P450 2B1 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from Rat Cytochrome P450 2B1 (301-400/491 aa)
Target	Cytochrome P450 2B1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid: 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	56 kDa
Uniprot ID	P00176
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human li...



Immunohistochemical analysis of formalin...



Immunohistochemical staining of human li...