


Product Datasheet

Cathepsin H heavy chain antibody (orb182888)

Description	Cathepsin H heavy chain antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Cathepsin H heavy chain (116-200/335 aa)
Target	Cathepsin H heavy chain
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	19 kDa
Uniprot ID	P09668
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse liver lys...