



Product Datasheet

DDB1 antibody (orb18292)

Description	Goat polyclonal antibody to DDB1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	DDB1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Concentration	0.01-0.03 µg/ml
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 140kDa band seen in HeLa, Jurkat and NSO (mouse) cell lysates. Recommended for use at 0.01-0.03 µg/ml.
Protein Sequence	DLIKVVEELTRIH
Clonality	Polyclonal
MW	127
Entrez	1642
Dilution Range	ELISA: 1:32000, WB: 0.01-0.03 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of NIH3T3 (A) and ...



Western blot analysis of HeLa (Lane 1), ...