



Product Datasheet

Cardiac Troponin T antibody (orb182936)


Description	Cardiac Troponin T antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Cardiac Troponin T (201-298/298 aa)
Target	Cardiac Troponin T
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	36 kDa
Uniprot ID	P45379
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, ICC=1:100-500, Flow-Cyt=1µg/Test, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat hear...



Immunohistochemical staining of human co...



Immunohistochemical staining of rat hear...