



Product Datasheet

FOXP2 antibody (orb18295)

Description	Goat polyclonal antibody to FOXP2
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	FOXP2 (C terminus)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000. WB: Approx 80kDa band observed in human brain (cerebellum) lysates (calculated MW of 79.9kDa according to NP_055306.1). Recommended concentration: 0.5-1.5 µg/ml.
Protein Sequence	REIEEEPLSEDLE
Clonality	Polyclonal
MW	79.9; 82.6; 70.1
Entrez	93986
Dilution Range	ELISA: 1:128000, WB: 0.5-2 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of human cerebellu...



Immunohistochemical staining of Human Te...