



Product Datasheet

MS4A1 antibody (orb18296)

Description	Goat polyclonal antibody to MS4A1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	CD20 / MS4A1 (C Terminus)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:4000.WB: Approx 35kDa band observed in Human Lymph lysates (calculated MW of 33.1kDa according to NP_068769.2). Recommended concentration: 0.003-0.01 µg/ml.
Protein Sequence	CQDQESSPIENDSSP
Clonality	Polyclonal
MW	33.1
Entrez	931
Dilution Range	ELISA: 1:128000, WB: 0.01-0.05 µg/ml, IHC-P: 10ug/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of human lymph lys...



Immunohistochemical staining of Human To...