

Product Datasheet

Cardiodilatin-related peptide antibody (orb182961)

Description	Cardiodilatin-related peptide antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human Cardiodilatin-related peptide (Processed active peptide) (26-55/153 aa)
Target	Cardiodilatin-related peptide
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	3/17 kDa
Uniprot ID	P01160
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.