



Product Datasheet

MAX antibody (orb18299)


Description	Goat polyclonal antibody to MAX
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Target	MAX
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 23kDa band seen in Jurkat Lysates. Recommended for use at 0.5-2 µg/ml.
Protein Sequence	EEPQSRKKLRMEAS
Clonality	Polyclonal
MW	18.3; 17.2
Entrez	4149
Dilution Range	ELISA: 1:32000, WB: 0.01-0.03 µg/ml, IHC-P: 10 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat lysate u...



Western blot analysis of HeLa lysate usi...



Immunohistochemical staining of Human Ad...