


Product Datasheet

NCF1 antibody (orb18301)


Description	Goat polyclonal antibody to NCF1
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, WB
Target	NCF1 / p47phox
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 45kDa band observed in lysates of cell line U937 (calculated MW of 44.7kDa according to NP_000256.3). Recommended concentration: 0.5-1.5 µg/ml.
Protein Sequence	SESTKRKLASAV
Clonality	Polyclonal
MW	44.7
Entrez	653361
Dilution Range	ELISA: 1:20000, WB: 0.01-1 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of U937 cell lysat...



Western blot analysis of Mouse Thymus (L...



Immunohistochemical staining of Human Co...