


Product Datasheet

MTR antibody (orb18307)

Description	Goat polyclonal antibody to MTR
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	MTR
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 150kDa band seen in HeLa and Jurkat cell lysates. Recommended for use at 1-3 µg/ml.
Protein Sequence	VEKWLGPILGYDTD
Clonality	Polyclonal
MW	141; 135. 96.1
Entrez	4548
Dilution Range	ELISA: 1:32000, IHC-P: 3-5 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Adipose Tissue (Adipocytes)