

Product Datasheet

DLC1 (phospho-Ser567) antibody (orb183119)

Description	DLC1 (phospho-Ser567) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human DLC1 around the phosphorylation site of Ser567 (PG(p-S)PD)
Target	phospho-DLC1 (Ser567)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	171 kDa
Uniprot ID	Q96QB1
Expiration Date	12 months from date of receipt.