


Product Datasheet

DLC1 antibody (orb183120)

Description	DLC1 antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Rho GTPase-activating protein 7 (1201-1300/1528 aa)
Target	DLC1 + DLC2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	171/122 kDa
Uniprot ID	Q96QB1
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...