



Product Datasheet

ITK antibody (orb18315)


Description	Goat polyclonal antibody to ITK
Species/Host	Goat
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Target	ITK
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 75 kDa bands seen in Jurkat Lysates. Recommended for use at 0.5-2 µg/ml.
Protein Sequence	RLLRQLAEIAESGL
Clonality	Polyclonal
MW	71.8
Entrez	3702
Dilution Range	ELISA: 1:2000, WB: 1-3g/ml, IHC-P: 10 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat nuclear ...



Immunofluorescence analysis of paraforma...



Immunofluorescence analysis of paraforma...