

Product Datasheet

DNAJC13 antibody (orb183156)

Description	DNAJC13 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated synthetic peptide derived from human DNAJC13 (11-120/2243 aa)
Target	DNAJC13
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	254 kDa
Uniprot ID	O75165
Expiration Date	12 months from date of receipt.