


Product Datasheet

FEM1C antibody (orb18319)

Description	Goat polyclonal antibody to FEM1C
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	FEM1C
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx 75kDa band observed in heart and kidney lysates from mouse and rat (calculated MW of 68.6kDa according to mouse NP_775599.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	HIPEKLETQVSLHR
Clonality	Polyclonal
MW	68.7
Entrez	56929
Dilution Range	ELISA: 1:8000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



orb18319
(0.5ug/ml)
staining of rat
kidn...