



Product Datasheet

PPP2R5E antibody (orb18322)

Description	Goat polyclonal antibody to PPP2R5E
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	PPP2R5E
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 55kDa and 40kDa bands seen in mouse brain extracts. Recommended for use at 0.5-2 µg/ml. Please note that the two bands may correspond to two different isoforms. Two isoforms of the epsilon mRNA have been detected in Northern blots (J Biol Chem 1996 Sep 6;271(36):22081-9).
Protein Sequence	LKRGLRRDGIPT
Clonality	Polyclonal
MW	54.7
Entrez	5529
Dilution Range	ELISA: 1:16000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



orb18322
staining
(0.5ug/ml) of
mouse br...



Western blot
analysis of
Human
Cerebellu...