



Product Datasheet

FAM3D antibody (orb183253)

Description	FAM3D antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human FAM3D (121-224/224 aa)
Target	FAM3D
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	22 kDa
Uniprot ID	Q96BQ1
Dilution Range	IHC-P=1:100-500, IHC-F=1:100-500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Mouse pl...



Immunohistochemical staining of rat sple...