

Product Datasheet

DILP5 antibody (orb183256)

Description	DILP5 antibody
Species/Host	Rabbit
Reactivity	Drosophila
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from Fire fly DILP5 (21-100/108 aa)
Target	DILP5
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	9.5 kDa
Uniprot ID	Q7KUD5
Expiration Date	12 months from date of receipt.