


Product Datasheet

SMAD4 antibody (orb18328)

Description	Goat polyclonal antibody to SMAD4
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IP
Target	SMAD4 / MADH4
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 60-65kDa band seen in mouse NIH-3T3 cell lysate. Recommended for use at 0.5-2 µg/ml.
Protein Sequence	HTMPIADPQPLD
Clonality	Polyclonal
MW	60.4
Entrez	4089
Dilution Range	ELISA: 1:32000
Expiration Date	12 months from date of receipt.



orb18328
staining
(0.5ug/ml) of
NIH-3T3 ...