


Product Datasheet

PPP2R5D antibody (orb18329)

Description	Goat polyclonal antibody to PPP2R5D
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, WB
Target	PPP2R5D
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:20000.WB: Approx 70kDa and 60kDa bands seen in 293 lysate (Predicted MW of 68kDa and 62kDa according to NP_006236 and AAH01175, respectively). Recommended for use at 1-3 µg/ml.
Protein Sequence	KRAEEFLTASQEAL
Clonality	Polyclonal
MW	70.0; 66.2; 58.5; 53.1
Entrez	5528
Dilution Range	ELISA: 1:20000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb18329
staining
(2ug/ml) of
293 lysate...