


Product Datasheet

PPP2R5B antibody (orb18330)

Description	Goat polyclonal antibody to PPP2R5B
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	B56 beta isoform
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:2000.WB: Approx. 55kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 57.4kDa according to NP_006235.1). Recommended concentration: 2-3 µg/ml.
Protein Sequence	QRLTPQVAASGGQS
Clonality	Polyclonal
MW	57.4
Entrez	5526
Dilution Range	ELISA: 1:2000, WB: 2-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Cerebellu...