

Product Datasheet

EEF1A2 (phospho-ser358) antibody (orb183311)

Description	EEF1A2 (phospho-ser358) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human EEF1A2 around the phosphorylation site of Ser358 (GY(p-S)PV)
Target	phospho-EEF1A2Ser358
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50 kDa
Uniprot ID	Q05639
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.