


Product Datasheet

ATF4 antibody (orb18333)

Description	Goat polyclonal antibody to ATF4
Species/Host	Goat
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	ATF4
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 40kDa band seen in Raji cell lysate. Recommended for use at 1-3 µg/ml.
Protein Sequence	EEVRKARGKKRVP
Clonality	Polyclonal
MW	38.6
Entrez	468
Dilution Range	ELISA: 1:128000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 lysate using ATF4 antibody.